ISSN 2046-1690



# Add Jaw-Thrust To Evolve 2-Rescuers Hands-Only CPR

Peer review status: No

*Corresponding Author:* Dr. Deepak Gupta, Anesthesiologist, Wayne State University, 48201 - United States of America

#### Submitting Author:

Dr. Deepak Gupta, Anesthesiologist, Wayne State University, 48201 - United States of America

*Other Authors:* Dr. Arvind Srirajakalidindi, Staff Anesthesiologist, Detroit Medical Center, Anesthesiology - United States of America

Article ID: WMC005593 Article Type: My opinion Submitted on:15-Sep-2019, 03:42:41 AM GMT Published on: 19-Sep-2019, 07:15:34 AM GMT Article URL: http://www.webmedcentral.com/article\_view/5593 Subject Categories:CARDIOLOGY Keywords:Jaw-Thrust, 2-Rescuers, Hands-Only CPR

**How to cite the article:**Gupta D, Srirajakalidindi A. Add Jaw-Thrust To Evolve 2-Rescuers Hands-Only CPR. WebmedCentral CARDIOLOGY 2019;10(9):WMC005593

**Copyright:** This is an open-access article distributed under the terms of the Creative Commons Attribution License(CC-BY), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Source(s) of Funding: NOT APPLICABLE

### Competing Interests:

NOT APPLICABLE

# Add Jaw-Thrust To Evolve 2-Rescuers Hands-Only CPR

Author(s): Gupta D, Srirajakalidindi A

### My opinion

Cardiocerebral resuscitation (CCR) as an evolved form of cardiopulmonary resuscitation (CPR) has been termed as Hands-Only CPR (HOCPR) for the lay rescuers.<sup>[1]</sup> iStan Adult Patient Simulator (CAE Healthcare, Saint-Laurent, Quebec, Canada) provides definitive evidence of "covert― pulmonary resuscitation happening during CCR. During a simulated clinical experience (SCE), end-tidal carbon dioxide (ETCO<sub>2</sub>) obtained at baseline was 25 mmHg (Figure 1A) when providing 788 ml tidal volume (TV) breaths during manual ventilation of intubated iStan. Subsequently obtained ETCO<sub>2</sub> was 2 mmHg (Figure 1B) when HOCPR generated 102 ml TV "breaths― in non-ventilated iStan with endotracheal tube open-to-room-air. Although delivery of prematurely early breaths during manual ventilation miscalculated respiratory rate (Figure 1A), minute ventilation (MV) equated to ~10 breaths/min during 7.3 I/min manual ventilation paradoxically increased to 9.9 I/min when ~100 compressions/min HOCPR generated 102 ml TV "breaths― (HOCPR-induced-high-frequency-"ventilation―). Essentially, the SCE with intubated iStan is making the case for keeping airway-and-mouth open during HOCPR. Therefore, HOCPR training should devise SCE for 2-rescuers HOCPR wherein second lay rescuers are providing jaw thrust with their elbows resting on the floor (Figure 1C-1D) to comfortably sustain jaw thrust and keeping airway-and-mouth open for prolonged periods during HOCPR. Moreover, future investigators should document if they observe that, after the end of each exhalation, there is evidence of measurable ETCO<sub>2</sub> as "cardiogenic oscillations― during each inspiratory downstroke in intubated and manually or mechanically ventilated patients during the times when they are receiving continuous chest compressions. Thereafter, they should explore whether controlling wind-speed around the open mouth can improve the "inhalability― of breathing zone room air in apneic mannequins.[2-3] Summarily, even though the above-mentioned observed evidence of pulmonary resuscitation during CCR (HOCPR) may need further exploration and validation studies, adding jaw thrust in the interim to

evolve 2-rescuers HOCPR may only contribute to bettering the process of CPR. $\hat{A}$ 

### References

- Ewy GA: Cardiocerebral and cardiopulmonary resuscitation - 2017 update. Acute Med Surg 2017; 4:227-34
- Schmees DK, Wu YH, Vincent JH: Visualization of the airflow around a life-sized, heated, breathing mannequin at ultralow windspeeds. *Ann Occup Hyg* 2008; 52:351-60
- Anthony TR, Anderson KR: Computational fluid dynamics investigation of human aspiration in low-velocity air: orientation effects on mouth-breathing simulations. *Ann Occup Hyg* 2013; 57:740-57

## Figure Legend

Figure 1: Evidence of Exhaled Carbon Dioxide (ETCO<sub>2</sub>) During HOCPR Prompting To Keep Airway And Mouth Open With Jaw Thrust During HOCPR [(1A): Baseline ETCO<sub>2</sub> at ~25 mmHg During Manual Ventilation; (1B): Evident ETCO<sub>2</sub> at ~2 mmHg During HOCPR; (1C): Schematically, Second Lay Rescuer Providing Jaw Thrust With Elbows Resting On Floor; (1D): Schematic Close Up View Demonstrating Jaw Thrust Providing Fingersâ€<sup>TM</sup> Alignment]

